


NEW MULTI-MYTHEN 3K DETECTOR EXCELLENT LaB_6 POWDER DATA COLLECTED IN JUST TWO MINUTES



Dr. Thomas Hartmann
Powder XRD and Sales
hartmann@stoe.com



A glass capillary of 0.5mm diameter was filled with NIST 660b LaB_6 standard material and mounted on a STOE STADI P in Debye-Scherrer geometry using a STOE MULTI MYTHEN2 R 3K detector (450 μm chip thickness) and Mo Ka1 radiation.

STOE MULTI-MYTHEN DETECTOR

MYTHEN2 R n K

(n = NUMBER OF MODULES WITH
 $n = 2$ OR 3)


STOE offers this new product for its STADI P and MP series to afford high resolution data collection in shortest measuring times

info@stoe.com.



Design study of the left- and right-handed STOE MYTHEN2 R 2K

STOE MYTHEN2 R 3K shown in picture above



The powder pattern above was measured from 2 to 110° 2 θ in 120 s, using the one step technique, and the data taken to refine the cell. The first 40° 2 θ of the fit are magnified in the upper right corner of the figure. The impressive results of the refinement are listed in the table below.

File title: LaB6
Wavelength: 0.709300
Number of accepted peaks: 68
2Theta window: 0.050
2Theta zeropoint: 0.0023 (refineable)
Symmetry: Cubic P
Spacegroup: P m-3m (No. 221)
Initial cell parameters: Cell_A: 4.1569

Refined cell parameters :
Cell_A: 4.15698(11)
Cell_Volume: 71.834(3)
Number of single indexed lines: 68
Number of unindexed lines: 0
2Theta zeropoint: 0.0014(23)
Average delta(2Theta) = 0.005
Figure of Merit F(30) = 460.4 (0.002, 31)